

Reporting Title: AMB/Candin Susceptibility Panel**Performing Location:** U.T. Health Science Ctr San Antonio**Shipping Instructions:**

1. See Infectious Specimen Shipping Guidelines.
2. Place specimen in a large infectious container (T146) and label as an etiologic agent/infectious substance.

Specimen Requirements:

Container/Tube: Agar slant or other appropriate media

Specimen Volume: Organism in pure culture

Complete and submit with specimen:

1. Fungus Testing Lab request form

Collection Instructions:

1. Organism must be in pure culture, actively growing.
2. Place specimen in a large infectious container (T146) and label as etiologic agent/infectious substance.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient		

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FAMCP	AMB/Candin Susceptibility Panel	Alphanumeric		Not Provided

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

Referral

CPT Code Information:

87999-Unlisted Microbiology Procedure (Refer to patient report to apply the appropriate CPT code in place of this unlisted CPT code based on specimen organism.)

87188 x 4 mould MIC microdilution or agar dilution (if appropriate)

87186 x 4 yeast MIC microdilution or agar dilution (if appropriate)