

# **Test Definition: FPPOS**

Posaconazole Susceptibility Testing

**Reporting Title:** Posaconazole (POS) Susceptibility **Performing Location:** U.T. Health Science Ctr San Antonio

# **Specimen Requirements:**

Container/Tube: Agar slant or other appropriate media

Specimen Volume: Organism in pure culture

Complete and submit with specimen:

1. Fungus Test Lab Request Form

#### Collection Instructions:

- 1. Organism must be in pure culture, actively growing.
- 2. Place specimen in a large infectious container (T146) and label as etiologic agent/infectious substance.

Specimen Type	Temperature	Time	Special Container	
Varies	Ambient			

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
FPPOS	Posaconazole (POS) Susceptibility	Alphanumeric		Not Provided

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

### **CPT Code Information:**

87999 - Unlisted Microbiology Procedure (Refer to patient report to apply the appropriate CPT code in place of this unlisted CPT code based on specimen organism.)

87188 - mould - MIC microdilution or agar dilution (if appropriate) 87186 - yeast - MIC microdilution or agar dilution (if appropriate)