

---

**Reporting Title:** Phencyclidine Confirmation, U  
**Performing Location:** Rochester**Ordering Guidance:**

1. For situations where chain of custody is required, a Chain of Custody Kit (T282) is available. For chain-of-custody testing, order PCPX / Phencyclidine Confirmation, Chain of Custody, Random, Urine
2. If urine creatinine is required or adulteration of the sample is suspected, order ADULT / Adulterants Survey, Random, Urine.
3. Additional drug panels and specific requests are available. Call 800-533-1710 or 507-266-5700.

**Specimen Requirements:**

Supplies: Urine Tubes, 10 mL (T068)

Container Tube: Plastic, 10 mL urine tube

Specimen Volume: 10 mL

Collection Instructions:

1. Collect a random urine specimen.
2. No preservative.

Additional Information:

1. No specimen substitutions.
2. STAT requests are not accepted for this test.

**Specimen Minimum Volume:**

1.2 mL

**Forms:**

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	72 hours	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
6673	Phencyclidine-by GC/MS	Alphanumeric	ng/mL	16254-5
21192	Phencyclidine Interpretation	Alphanumeric		69050-3
21193	Chain of Custody	Alphanumeric		77202-0

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

83992  
G0480 (if appropriate)

**Reference Values:**

Negative  
Cutoff concentrations  
Phencyclidine by GC-MS: <10 ng/mL