

Pipecolic Acid, Serum

# **Reporting Title:** Pipecolic Acid, S **Performing Location:** Rochester

### **Necessary Information:**

Patient's age is required.

#### **Specimen Requirements:**

Patient Preparation: Fasting 12 hours or more. (Collect specimens from infants and small children just before next feeding) Supplies: Sarstedt Aliquot Tube, 5 mL (T914) Collection Container/Tube: Preferred: Serum gel Acceptable: Red top Submission Container/Tube: Plastic vial Specimen Volume: 1.5 mL Collection Instructions: Centrifuge and aliquot serum into plastic vial.

#### **Specimen Minimum Volume:**

1 mL

### Forms:

If not ordering electronically, complete, print, and send a Biochemical Genetics Test Request (T798) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	94 days	
	Refrigerated	14 days	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
81326	Pipecolic Acid, S	Numeric	nmol/mL	32334-5
29962	Interpretation	Alphanumeric		59462-2
29964	Reviewed By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

## Supplemental Report:



Pipecolic Acid, Serum

No

# **CPT Code Information:**

82542

## **Reference Values:**

<6 months: < or =6.0 nmol/mL 6 months-<1 year: < or =5.9 nmol/mL 1-17 years: < or =4.3 nmol/mL > or =18 years: < or =7.4 nmol/mL