

# **Test Definition: UBT**

Helicobacter pylori Breath Test

Reporting Title: H. pylori C Urea Breath Test

Performing Location: Rochester

### **Ordering Guidance:**

An alternative test for the diagnosis of active Helicobacter pylori infection in patients is the HPFRP / Helicobacter pylori with Clarithromycin Resistance Prediction, Molecular Detection, PCR, Feces, which requires a different collection.

### **Specimen Requirements:**

Patient Preparation:

- 1. Patient should be fasting for 1 hour.
- 2. Do not administer this test if this list of instructions is not followed, as test results may be inaccurate:
- a. Patients should not take proton-pump inhibitors (eg, Prilosec, Prevacid, Aciphex, Protonix, and Nexium) or bismuth compounds (eg, Pepto-Bismol) for 2 weeks prior to testing.
- b. Patients should not take antibiotics for 4 weeks prior to testing.
- c. Histamine 2-receptor antagonists (H2RA), such as Pepcid, Tagamet, Axid, or Zantac may impact the sensitivity of the test and, if possible, should be discontinued for 2 weeks prior to testing.
- 3. Carafate (sucralfate) does not interfere with the test. Use of antacids does not affect the accuracy of this assay.

Supplies: H. Pylori Breath Kit - Meridian Breath ID (T906: fees apply)

Collection Instructions:

- 1. Do not collect if patient is younger than 3 years of age.
- 2. Follow instructions included with kit.
- 3. Mixing the (13)C-Urea Tablet
- a. Dissolve the Citrica and the (13)C-enriched urea tablet in 150 to 200 mL (5.1 to 6.8 oz.) of tap water in the provided drinking cup.
- b. Close the lid firmly using both hands. Place fingers over lid and shake thoroughly for a few minutes, until the Citrica Powder and the urea tablet are completely dissolved.

Note: Tiny particles may remain visible after thorough mixing. However, if more substantial particulate matter is still present after five minutes of mixing, discard the solution and repeat the procedure with a new kit.

## **Specimen Minimum Volume:**

Bag of "breath" must be full

#### Forms:

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

- -Infectious Disease Serology Test Request (T916)
- -Gastroenterology and Hepatology Test Request (T728)
- -Microbiology Test Request (T732)
- -General Request (T239)



# **Test Definition: UBT**

Helicobacter pylori Breath Test

Specimen Type	Temperature	Time	Special Container	
Breath	Ambient	14 days	BREATH TEST BAG	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
81590	H. pylori C Urea Breath Test	Alphanumeric		29891-9

LOINC and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

### **CPT Code Information:**

83013

## **Reference Values:**

Negative

Reference values apply to all ages.