**Test Definition: GRAS2** 

Grass Panel # 2, Serum

**Reporting Title:** Grass Panel # 2 **Performing Location:** Rochester

### Ordering Guidance:

This test uses a pooled allergen reagent; therefore, the multi-allergen Immunocap (panel cap) is reported with a single qualitative class result and concentration. This is the appropriate first-tier test for allergic disease.

For a listing of allergens available for testing, see Allergens - Immunoglobulin E (IgE) Antibodies

# **Specimen Requirements:**

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL for every 5 allergens requested

Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

### **Specimen Minimum Volume:**

0.3 mL

#### Forms:

If not ordering electronically, complete, print, and send an Allergen Test Request (T236) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	90 days	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
GRAS2	Grass Panel # 2	Numeric	kU/L	30188-7

LOINC and CPT codes are provided by the performing laboratory.

### Supplemental Report:

No



**Test Definition: GRAS2** 

Grass Panel # 2, Serum

## **CPT Code Information:**

86003

## **Reference Values:**

Class	IgE kU/L	Interpretation
0	<0.10	Negative
0/1	0.10-0.34	Borderline/equivocal
1	0.35-0.69	Equivocal
2	0.70-3.49	Positive
3	3.50-17.4	Positive
4	17.5-49.9	Strongly positive
5	50.0-99.9	Strongly positive
6	> or =100	Strongly positive

Reference values apply to all ages.