

**Reporting Title:** Thyroid Autoantibodies Profile, S  
**Performing Location:** Rochester

**Ordering Guidance:**

If thyroglobulin tumor marker testing is desired, order HTG2 / Thyroglobulin, Tumor Marker, Serum.

**Specimen Requirements:**

Patient Preparation: For 12 hours before specimen collection do not take multivitamins or dietary supplements containing biotin (vitamin B7), which is commonly found in hair, skin, and nail supplements and multivitamins.

Collection Container/Tube: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

**Specimen Minimum Volume:**

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Frozen	30 days	
	Ambient	7 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
TPO	Thyroperoxidase Ab, S <b>Also used by tests: TPO</b>	Numeric	IU/mL	8099-4
TGAB	Thyroglobulin Antibody, S <b>Also used by tests: TGAB</b>	Numeric	IU/mL	56536-6

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

---

**Components:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
TPO	Thyroperoxidase Ab, S			Yes	Yes
TGAB	Thyroglobulin Antibody, S			Yes	Yes

**CPT Code Information:**

86376-Thyroperoxidase antibody

86800-Thyroglobulin antibody

**Reference Values:**

## THYROGLOBULIN ANTIBODY

&lt;4.0 IU/mL

Reference values apply to all ages.

## THYROPEROXIDASE ANTIBODIES

&lt;9.0 IU/mL

Reference values apply to all ages.