

# **Test Definition: SCOC**

Coccidioides Antibody, Complement Fixation and Immunodiffusion, Serum

# **Reporting Title:** Coccidioides Ab, CompF/ImmDiff, S **Performing Location:** Rochester

#### Ordering Guidance:

This test is a confirmatory assay for positive screening tests (ie, enzyme immunoassay).

The recommended test to evaluate for possible Coccidioides infection is COXIS / Coccidioides Antibody Screen with Reflex, Serum.

#### Specimen Requirements:

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top Submission Container/Tube: Plastic vial Specimen Volume: 1.8 mL Collection Instructions: Centrifuge and aliquot serum into plastic vial.

#### **Specimen Minimum Volume:**

1.2 mL

#### Forms:

If not ordering electronically, complete, print, and send Infectious Disease Serology Test Request (T916) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	14 days	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
8295	Coccidioides Ab, CompF, S	Alphanumeric		In Process
21649	Coccidioides, IgG, ImmDiff, S	Alphanumeric		62459-3
21648	Coccidioides, IgM, ImmDiff, S	Alphanumeric		62458-5

LOINC and CPT codes are provided by the performing laboratory.



# **Test Definition: SCOC**

Coccidioides Antibody, Complement Fixation and Immunodiffusion, Serum

# Supplemental Report:

No

## **CPT Code Information:**

86635 x 3

### **Reference Values:**

COMPLEMENT FIXATION Negative If positive, results are titered.

IMMUNODIFFUSION Negative Results are reported as positive, negative, or equivocal.