

# **Test Definition: HEAG**

Hepatitis B Virus e Antigen and Hepatitis B Virus e Antibody, Serum

Reporting Title: Hepatitis Be Ag and Ab, S

Performing Location: Rochester

### **Ordering Guidance:**

If ordered with HBVQN / Hepatitis B Virus (HBV) DNA Detection and Quantification by Real-Time PCR, Serum; send separate vials.

## **Necessary Information:**

Date of collection is required

## Specimen Requirements:

Patient Preparation: For 24 hours before specimen collection, patient should not take multivitamins or dietary supplements (eg, hair, skin, and nail supplements) containing biotin (vitamin B7).

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Serum gel Submission Container/Tube: Plastic vial

Specimen Volume: 0.8 mL Collection Instructions:

1. Centrifuge blood collection tube per manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).

2. Aliquot serum into a plastic vial.

### Specimen Minimum Volume:

0.7 mL

#### Forms:

If not ordering electronically, complete, print, and send 1 of the following:

- -Gastroenterology and Hepatology Test Request (T728)
- -Infectious Disease Serology Test Request (T916)

Specimen Type	Temperature	Time	Special Container
Serum SST	Frozen (preferred)	90 days	
	Refrigerated	6 days	
	Ambient	72 hours	



# **Test Definition: HEAG**

Hepatitis B Virus e Antigen and Hepatitis B Virus e Antibody, Serum

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
EAG	Hepatitis Be Ag, S	Alphanumeric		13954-3
	Also used by tests: EAG			
HEAB	HBe Antibody, S	Alphanumeric		33463-1
	Also used by tests: HEAB			

LOINC and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

# **Components:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
EAG	Hepatitis Be Ag, S			Yes	Yes
HEAB	HBe Antibody, S			Yes	Yes

### **CPT Code Information:**

86707 87350

### **Reference Values:**

HEPATITIS Be ANTIGEN: Negative

**HEPATITIS Be ANTIBODY:** 

Negative

For more information see Viral Hepatitis Serologic Profiles.