

**Reporting Title:** Drug Screen, Prescription/OTC, S**Performing Location:** Rochester**Ordering Guidance:**

This test is not performed using chain of custody. For chain-of-custody testing, order DSSX / Drug Screen, Prescription/Over the Counter, Chain of Custody, Serum.

**Specimen Requirements:**

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Red top (serum gel/SST tubes are not acceptable)

Submission Container/Tube: Plastic vial

Specimen Volume: 2.75 mL

Collection Instructions: Centrifuge and aliquot serum into plastic vial within 2 hours of collection.

**Specimen Minimum Volume:**

1.1 mL

**Forms:**

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	3 hours	

**Ask at Order Entry (AOE) Questions:**

Test ID	Question ID	Description	Type	Reportable
DSS	SU_DR	Suspect Drug	Plain Text	No

---

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
31072	Drugs detected:	Alphanumeric		20785-2
31168	Chain of Custody	Alphanumeric		77202-0

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**  
No**CPT Code Information:**

80307

**Reference Values:**

Drugs detected are presumptive. Additional testing may be required to confirm the presence of any drugs detected.