

Reporting Title: Flurazepam (Dalmane)**Performing Location:** Medtox Laboratories, Inc.**Specimen Requirements:**

Submit only 1 of the following specimens:

Plasma

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 2 mL sodium heparin plasma frozen in plastic vial.

Serum

Draw blood in a plain, red-top tube(s), serum gel tube is not acceptable. Spin down and send 2 mL of serum frozen in plastic vial.

Specimen Minimum Volume:

0.3 mL

Specimen Type	Temperature	Time	Special Container
Varies	Frozen	180 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z1133	Flurazepam	Alphanumeric		46979-1
Z1160	Desalkylflurazepam	Alphanumeric		46980-9
Z1161	Flurazepam + Desalkylflurazepam	Alphanumeric		73826-0

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80346

Reference Values:

Flurazepam:

Reference Range: 0 - 30 ng/mL

Desalkylflurazepam:

Reference Range: 30 - 150 ng/mL

Flurazepam + Desalkylflurazepam:

Reference Range: 30 - 180 ng/mL