## **Test Definition: FDM**

Dextromethorphan (DM), Serum

**Reporting Title:** Dextromethorphan (DM) **Performing Location:** Medtox Laboratories, Inc.

## **Specimen Requirements:**

Submit only 1 of the following specimens:

#### Plasma

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 2 mL sodium heparin refrigerated in a plastic vial.

#### Serum

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

## **Specimen Minimum Volume:**

0.25 mL

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z1148	Dextromethorphan (DM)	Alphanumeric		13623-4

LOINC and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

#### **CPT Code Information:**

80362



# **Test Definition: FDM**

Dextromethorphan (DM), Serum

#### **Reference Values:**

Reference Range: 2.0 - 6.0 ng/mL