

Vigabatrin (Sabril)

## **Reporting Title:** Vigabatrin (Sabril)

Performing Location: Medtox Laboratories, Inc.

#### **Specimen Requirements:**

Submit only 1 of the following specimens

Plasma

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 1 mL sodium heparin plasma refrigerated in a plastic vial.

Serum

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

### **Specimen Minimum Volume:**

0.25 mL

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	
	Ambient	72 hours	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z1176	Vigabatrin (Sabril)	Alphanumeric		30042-6

LOINC and CPT codes are provided by the performing laboratory.

## Supplemental Report:

No

## **CPT Code Information:**

80339



# **Test Definition: VIGA**

Vigabatrin (Sabril)

**Reference Values:**