

Diltiazem (Cardizem, Dilacor)

# Reporting Title: Diltiazem (Cardizem, Dilacor)

Performing Location: Medtox Laboratories, Inc.

#### **Specimen Requirements:**

Submit only 1 of the following specimens:

Plasma

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 3 mL sodium heparin plasma refrigerated in a plastic vial.

Serum

Draw blood in a plain, red-top tube(s), serum gel tube is not acceptable. Spin down and send 3 mL of serum refrigerated in a plastic vial.

#### **Specimen Minimum Volume:**

0.6 mL

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	
	Ambient	72 hours	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z1188	Diltiazem	Alphanumeric		6898-1

LOINC and CPT codes are provided by the performing laboratory.

# Supplemental Report:

No

## **CPT Code Information:**

80299



Diltiazem (Cardizem, Dilacor)

#### **Reference Values:**

Reference Range: 50 - 200 ng/mL