

## **Test Definition: FAMP**

Amphetamine, Serum or Plasma

Reporting Title: Amphetamine, Serum or Plasma

**Performing Location:** Medtox Laboratories, Inc.

## **Specimen Requirements:**

Submit only 1 of the following specimens:

#### Plasma

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 5 mL sodium heparin plasma refrigerated in a plastic vial.

#### Serum

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 5 mL of serum refrigerated in a plastic vial.

### **Specimen Minimum Volume:**

0.6 mL

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Varies        | Refrigerated (preferred) | 14 days  |                   |
|               | Frozen                   | 180 days |                   |

#### **Result Codes:**

| Result ID | Reporting Name | Туре         | Unit | LOINC®  |
|-----------|----------------|--------------|------|---------|
| Z1206     | Amphetamine    | Alphanumeric |      | 42576-9 |

LOINC and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

#### **CPT Code Information:**

80324



# **Test Definition: FAMP**

Amphetamine, Serum or Plasma

**Reference Values:**