

## **Test Definition: FAVI**

Hypersensitivity Pneumonitis Avian Panel

Reporting Title: Avian Panel

Performing Location: Medical College of WI Allergy-Immunology MACC Fund Research Ctr.

## **Specimen Requirements:**

Collection Container/Tube: Red Top or SST Submission Container/Tube: Plastic vial

Specimen Volume: 3 mL Collection Instructions:

1. Draw blood in a Red-top tube(s). Serum-gel tube(s) is also acceptable.

2. Centrifuge and send 3 mL of serum refrigerated.

## **Specimen Minimum Volume:**

1 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	30 days	
	Frozen	365 days	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z0288	Pigeon Sera	Alphanumeric		Not Provided
Z0289	Pigeon DE	Alphanumeric		Not Provided
Z0286	Cockatiel	Alphanumeric		Not Provided
Z0285	Parakeet	Alphanumeric		Not Provided
Z0287	Parrot	Alphanumeric		Not Provided

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

## **CPT Code Information:**

86331x5



# **Test Definition: FAVI**

Hypersensitivity Pneumonitis Avian Panel

## **Reference Values:**

This panel includes the following antigens:

Pigeon Sera Pigeon DE Cockatiel Parakeet Parrot

Results are reported as: Negative Weakly Positive Positive

Strong Positive

This result must be correlated with patient's clinical response and should not solely be considered in the diagnosis.