

**Reporting Title:** Pepsinogen I**Performing Location:** Inter Science Institute**Specimen Requirements:**

Patient preparation:

Patient should be fasting 10-12 hours prior to collection of specimen. Antacids or other medications affecting stomach acidity or gastrointestinal motility should be discontinued, if possible, for at least 48 hours prior to collection.

Specimen Type: Serum

Container/Tube: Red top or SST

Specimen Volume: 3 mL

Collection Instructions: Draw blood in a plain, red-top tube(s), serum-gel tube(s) is acceptable. Separate immediately and send 3 mL of serum frozen in a plastic vial.

**Specimen Minimum Volume:**

1 mL

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	14 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
Z0941	Pepsinogen I	Alphanumeric	ng/mL	2736-7

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

83520

---

**Reference Values:**

24 - 214 ng/mL