

Reporting Title: Qualitative Fructose,Semen**Performing Location:** Rochester**Specimen Requirements:**

Patient Preparation: Patient should have 2 to 7 days of sexual abstinence at the time of semen collection.

Submit only 1 of the following specimens:

Specimen Type: Semen

Collection Container/Tube: Sterile container

Submission Container/Tube: Plastic container

Specimen Volume: Total ejaculate

Collection Instructions: Do not dilute specimen. Freeze specimen at -20 degrees C.

Specimen Type: Seminal plasma

Container/Tube: Plastic vial

Specimen Volume: 0.5 mL

Collection Instructions:

1. After semen collection, wait 30 to 40 minutes until the semen is liquefied, then centrifuge the semen for 10 minutes at maximum centrifuge speed.
2. Remove top 3/4 of specimen with a pipet and place in a plastic vial. Freeze specimen at -20 degrees C.
3. Discard remainder of centrifuged specimen.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container
Semen	Frozen		

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FRUCS	Fructose, Semen/Seminal P	Alphanumeric		13943-6
CLSM	Collection Site	Alphanumeric		56816-2
VOLSM	Semen Volume	Numeric	mL	3160-9
CMT49	Comment	Alphanumeric		48767-8

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82757

Reference Values:

Positive