

**Reporting Title:** Cytoplasmic Neutrophilic Ab, S**Performing Location:** Rochester**Additional Testing Requirements:**

When used for diagnosis, it is recommended that specific tests for proteinase 3, antineutrophil cytoplasmic antibodies (ANCA), and myeloperoxidase ANCA be performed in addition to testing for cytoplasmic ANCA and perinuclear ANCA.(1) This panel of tests is available; order VASC / Antineutrophil Cytoplasmic Antibodies Vasculitis Panel, Serum.

**Specimen Requirements:**

Collection Container/Tube:

Preferred: Serum gel

Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 0.8 mL

Collection Instructions: Centrifuge and aliquot serum into plastic vial.

**Specimen Minimum Volume:**

0.4 mL

**Forms:**

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

-General Request (T239)

-Renal Diagnostics Test Request (T830)

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	21 days	
	Frozen	21 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
3114	c-ANCA	Alphanumeric		In Process
3119	p-ANCA	Alphanumeric		17357-5

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

---

No

**CPT Code Information:**

86036 x2

86037-Titer (if appropriate)

**Reference Values:**

Negative

If positive for antineutrophil cytoplasmic antibodies, results are titered.