

Test Definition: APTSC

Activated Partial Thromboplastin Time (APTT),

Plasma

Reporting Title: Activated Partial Thrombopl Time, P

Performing Location: Rochester

Necessary Information:

Heparin or Coumadin therapy should be noted.

Specimen Requirements:

Only orderable as part of a special coagulation profile or as a reflex. For more information see:

ALUPP / Lupus Anticoagulant Profile, Plasma

ALBLD / Bleeding Diathesis Profile, Limited, Plasma

AATHR / Thrombophilia Profile, Plasma and Whole Blood

APROL / Prolonged Clot Time Profile, Plasma

ADIC / Disseminated Intravascular Coagulation/Intravascular Coagulation and Fibrinolysis (DIC/ICF) Profile, Plasma

Specimen Type	Temperature	Time	Special Container	
Plasma Na Cit	Frozen	14 days		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
APTSC	Activated Partial Thrombopl Time, P	Numeric	sec	14979-9

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

85730

Reference Values:

Only orderable as part of a special coagulation profile or as a reflex. For more information see:

ALUPP / Lupus Anticoagulant Profile, Plasma

ALBLD / Bleeding Diathesis Profile, Limited, Plasma

AATHR / Thrombophilia Profile, Plasma and Whole Blood

APROL / Prolonged Clot Time Profile, Plasma

ADIC / Disseminated Intravascular Coagulation/Intravascular Coagulation and Fibrinolysis (DIC/ICF) Profile, Plasma

25-37 seconds

The APTT may be 35% longer in full-term newborns that reach adult reference range by age 3 months and twice the adult upper limit in premature infants reaching adult reference range by age 6 months.