

**Test Definition: 22C3** 

Programmed Death-Ligand 1 (PD-L1) (22C3), Semi-Quantitative Immunohistochemistry,

Manual

Reporting Title: PD-L1 (22C3) SemiQuant IHC, Manual

Performing Location: Rochester

### **Ordering Guidance:**

For information on selection of programmed cell death 1-ligand 1 (PD-L1) testing, see <a href="https://news.mayocliniclabs.com/pd-l1-by-immunohistochemistry/">https://news.mayocliniclabs.com/pd-l1-by-immunohistochemistry/</a>.

## **Shipping Instructions:**

Attach the green pathology address label included in the kit to the outside of the transport container.

#### **Necessary Information:**

A pathology/diagnostic report and a brief history, including primary site of neoplasm, are required.

# **Specimen Requirements:**

Specimen Type: Tissue

Supplies: Pathology Packaging Kit (T554)

Collection Instructions: Formalin-fixed, paraffin-embedded tissue block; or 3 unstained glass, "positively charged" slides

with 4-microns formalin-fixed, paraffin-embedded tissue

Additional Information: One slide will be stained with hematoxylin and eosin and returned.

#### Forms:

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen: Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains Request (T763)

Oncology Test Request (T729)

Specimen Type	Temperature	Time	Special Container
Special	Ambient (preferred)		
	Refrigerated		

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
603763	Interpretation	Alphanumeric		59465-5
603764	Participated in the Interpretation	Alphanumeric		No LOINC Needed
603765	Report electronically signed by	Alphanumeric		19139-5
603766	Material Received	Alphanumeric		81178-6
603767	Disclaimer	Alphanumeric		62364-5
603768	Case Number	Alphanumeric		80398-1

LOINC® and CPT codes are provided by the performing laboratory.



# **Test Definition: 22C3**

Programmed Death-Ligand 1 (PD-L1) (22C3), Semi-Quantitative Immunohistochemistry, Manual

Supplemental Rep	port:	:
------------------	-------	---

No

## **CPT Code Information:**

88360