

Test Definition: HK27M

Histone H3 K27M Mutant (H3 K27M) Immunostain, Technical Component Only

Reporting Title: Histone H3 K27M IHC, Tech Only

Performing Location: Rochester

Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

Mayo Clinic Laboratories has multiple histone immunostains available. See table for ordering guidance.

			Mayo Clinic
Test ID	Published name	Indication	slide label
HG34	Histone 3.3 G34W (H3F3A G34W)	Giant cell tumor	H3 G34W
w	Immunostain, Technical Component Only	of bone (GCTB)	
HK27M	Histone H3 K27M Mutant (H3 K27M)	K27M mutant	H3 K27M
	Immunostain, Technical Component Only	midline gliomas	
HISME	Histone H3 Trimethyl K27 (H3 K27me[3])	MPNST and	H3 K27me3
	Immunostain, Technical Component Only	K27M mutant	
		midline gliomas	
НК36М	Histone H3 K36M Mutant (H3F3 K36M)	Chondroblastom	H3 K36M
	Immunostain, Technical Component Only	a	

Shipping Instructions:

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

Specimen Requirements:

Supplies: Immunostain Technical Only Envelope (T693)

Specimen Type: Tissue

Container/Tube: Immunostain Technical Only Envelope

Preferred: 2 Unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.

Acceptable: Formalin-fixed, paraffin-embedded tissue block

Forms:

If not ordering electronically, complete, print, and send an <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains</u>
<u>Request</u> (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
604990	Histone H3 K27M IHC, Tech Only	Alphanumeric		Bill only; no result



Test Definition: HK27M

Histone H3 K27M Mutant (H3 K27M) Immunostain, Technical Component Only

 $\ensuremath{\mathsf{LOINC}}\xspace^{\ensuremath{\texttt{@}}}$ and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

88342-TC, primary 88341-TC, if additional IHC

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
IHTOI	IHC Initial, Tech Only	1	88342	No	No
IHTOA	IHC Additional, Tech Only	1	88341	No	No