

Test Definition: UHSD2

Hemosiderin, Random, Urine

Reporting Title: Hemosiderin, Random, U

Performing Location: Rochester

Specimen Requirements:

Container/Tube: Plastic urine container

Specimen Volume: 13 mL **Collection Instructions:**

1. Collect a random urine specimen.

2. No preservative.

Specimen Type	Temperature	Time	Special Container	
Urine	Refrigerated (preferred)	7 days		
	Frozen	7 days		
	Ambient	2 hours		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
HSDU2	Hemosiderin, Random, U	Alphanumeric		4644-1

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83070

Reference Values:

Negative