

## **Test Definition: CARNP**

Carbapenemase Detection-Carba NP Test, Varies

# **Reporting Title:** Carbapenemase-Carba NP Test **Performing Location:** Rochester

#### Shipping Instructions:

- 1. For shipping information see Infectious Specimen Shipping Guidelines.
- 2. Place specimen in a large infectious container (T146) and label as an etiologic agent/infectious substance.

#### **Necessary Information:**

Specimen source and organism identification are required.

#### **Specimen Requirements:**

Specimen Type: Organism
Supplies: Infectious Container, Large (T146)
Container/Tube: Slant
Specimen Volume: Isolate
Collection Instructions: Submit Enterobacterales or *Pseudomonas aeruginosa* isolate in pure culture (ie, not mixed with other organisms), actively growing.

#### Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T244) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Refrigerated		

#### Ask at Order Entry (AOE) Questions:

Test	Question	Description	Туре	Reportabl
ID	ID			е
CARNP	Q00M0076	Specimen Source (Required) and Organism Identification	Plain	Yes
		(Required)	Text	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
CARNP	Carbapenemase-Carba NP Test	Alphanumeric	mL	74676-8

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No



### **Test Definition: CARNP**

Carbapenemase Detection-Carba NP Test, Varies

#### **Reflex Tests:**

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
CARNB	Carbapenemase-Carba NP Test	1	87185	Yes	No

#### **Reference Values:**

Negative