

# **Test Definition: LECT2**

Leukocyte Cell-Derived Chemotaxin 2 (LECT2), Immunostains Without Interpretation

Reporting Title: LECT2 IHC, Tech Only

Performing Location: Rochester

#### **Ordering Guidance:**

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

#### **Shipping Instructions:**

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

## **Specimen Requirements:**

Supplies: Immunostain Technical Only Envelope (T693)

**Specimen Type:** Tissue

Container/Tube: Immunostains Technical Only Envelope

**Preferred:** 2 Unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.

Acceptable: Formalin-fixed, paraffin-embedded tissue block

#### Forms:

If not ordering electronically, complete, print, and send a <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains</u>
Request (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
70809	LECT2 IHC, Tech Only	Alphanumeric		Bill only; no result

LOINC® and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

#### **CPT Code Information:**

88342-TC, primary 88341-TC, if additional IHC

#### **Reflex Tests:**

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
---------	----------------	-----------	----------	------------------	----------------------



# **Test Definition: LECT2**

Leukocyte Cell-Derived Chemotaxin 2 (LECT2), Immunostains Without Interpretation

IHTOI	IHC Initial, Tech Only	1	88342	No	No
IHTOA	IHC Additional, Tech Only	1	88341	No	No

# **Reference Values:**