

Test Definition: FADFL

Adenosine Deaminase in Peritoneal Fluid

Reporting Title: Adenosine Deaminase Peritoneal Fld

Performing Location: ARUP Laboratories

Specimen Requirements:

Specimen Type: Peritoneal fluid (Ascites, Paracentesis)

Container/Tube: Standard transport container

Specimen volume: 0.5 mL

Collection Instructions: Collect Peritoneal Fluid in a leak-proof container. Centrifuge specimen at room temperature, transfer 0.5 mL peritoneal fluid to plastic vial and Ship frozen. The specimen must remain frozen until received at the

performing lab.

Specimen Type	Temperature	Time	Special Container	
Peritoneal	Frozen (preferred)	30 days		
	Refrigerated	7 days		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
FADFL	Adenosine Deaminase Peritoneal Fld	Alphanumeric		49759-4

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84311

Reference Values:

0 - 30 U/L