

Test Definition: FADCF

Adenosine Deaminase in CSF

Reporting Title: Adenosine Deaminase in CSF

Performing Location: ARUP Laboratories

Specimen Requirements:

Collect CSF in a leak-proof container. Centrifuge specimen at room temperature and send 0.5 mL frozen.

Specimen Type	Temperature	Time	Special Container	
CSF	Frozen (preferred)	30 days		
	Refrigerated	7 days		
	Ambient	24 hours		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z5892	Adenosine Deaminase, CSF	Alphanumeric		35703-8

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84311

Reference Values:

0-9 U/L