

Test Definition: PRA

Renin Activity, Plasma

Reporting Title: Renin Activity, P
Performing Location: Rochester

Specimen Requirements:

Patient Preparation: For 4 to 6 weeks before specimen collection, spironolactone (Aldactone) should be discontinued, as plasma renin activity cannot be interpreted if the patient is being treated with spironolactone.

Collection Container/Tube: Chilled, lavender top (EDTA)

Submission Container/Tube: Plastic vial

Specimen Volume: 0.8 mL Collection Instructions:

- 1. Draw blood in a chilled syringe from a patient in a seated position; place specimen in a chilled, 3-mL lavender-top (EDTA) tube and mix well. Alternatively, draw blood directly in a chilled, EDTA tube.
- 2. Immediately place EDTA tube into an ice-water bath until thoroughly cooled.
- 3. Immediately centrifuge using a refrigerated centrifuge and aliquot plasma into a plastic vial.

Note: If a refrigerated centrifuge is unavailable, chill the centrifuge carriers before use. Centrifuge specimen for 5 minutes or less, then promptly aliquot plasma into a plastic vial.

4. Immediately freeze plasma.

Forms:

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

- -General Request (T239)
- -Cardiovascular Test Request (T724)

| Specimen Type | Temperature | Time | Special Container |
|---------------|-------------|---------|-------------------|
| Plasma EDTA | Frozen | 14 days | |

Result Codes:

| Result II |) | Reporting Name | Туре | Unit | LOINC® |
|-----------|---|-------------------|---------|---------|--------|
| 8060 | | Renin Activity, P | Numeric | ng/mL/h | 2915-7 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84244

Reference Values:

0-2 years: 4.6 ng/mL/h (mean)* Range: 1.4-7.8 ng/mL/h 3-5 years: 2.5 ng/mL/h (mean)* Range: 1.5-3.5 ng/mL/h 6-8 years: 1.4 ng/mL/h (mean)* Range: 0.8-2.0 ng/mL/h 9-11 years: 1.9 ng/mL/h (mean)* Range: 0.9-2.9 ng/mL/h



Test Definition: PRA

Renin Activity, Plasma

12-17 years: 1.8 ng/mL/h (mean)* Range: 1.2-2.4 ng/mL/h

Mean data not standardized as to time of day or diet. Infants were supine, children sitting.

*Stalker HP, Holland NH, Kotchen JM, Kotchen TA. Plasma renin activity in healthy children. J Pediatr.

1976;89(2):256-258

Na-depleted, upright (peripheral vein specimen)

18-39 years: 10.8 ng/mL/h (mean)

2.9-24.0 ng/mL/h (range)

> or =40 years: 5.9 ng/mL/h (mean)

2.9-10.8 ng/mL/h (range)

Na-replete, upright (peripheral vein specimen)

18-39 years: 1.9 ng/mL/h (mean) < or =0.6-4.3 ng/mL/hour (range) > or =40 years: 1.0 ng/mL/h (mean) < or =0.6-3.0 ng/mL/h (range)