

Lead, Nails

# Reporting Title: Lead, Nails

Performing Location: Rochester

# **Necessary Information:**

Indicate source of nails (fingernails or toenails), if known

# **Specimen Requirements:**

Supplies: Hair and Nails Collection Kit (T565)

Specimen Volume: 0.2 g

#### **Collection Instructions:**

- 1. Prepare and transport specimen per the instructions in the kit or see <u>Collecting Hair and Nails for Metals Testing</u>.
- 2. Clippings should be taken from all 10 fingernails or toenails.

Specimen Type	Temperature	Time	Special Container
Nail	Ambient (preferred)		
	Frozen		
	Refrigerated		

# Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
PBNA	PBNSC	Specimen Source:	Answer List	Yes
		• Finger		
		• Toe		
		<ul> <li>Finger and Toe</li> </ul>		

# **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
2506	Lead, Nails	Numeric	mcg/g	8202-4
PBNSC	Specimen Source	Alphanumeric		31208-2

 $\ensuremath{\mathsf{LOINC}}\xspace^{\ensuremath{\mathsf{\$}}}$  and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

**CPT Code Information:** 83655

Reference Values:

<4.0 mcg/g of nails Reference values apply to all ages.



# **Test Definition: PBNA**

Lead, Nails