

Test Definition: PINW

Pinworm Exam, Perianal

Reporting Title: Pinworm Exam, Perianal

Performing Location: Rochester

Specimen Requirements:

Supplies: Swubes (T300) **Specimen Type:** Perianal

Container/Tube: SWUBE disposable paddle (Falcon) or similar method of collection

Specimen Volume: Entire specimen

Collection Instructions: See <u>Pinworm Collection Instructions</u>

Forms:

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

- -Microbiology Test Request (T244)
- -Gastroenterology and Hepatology Test Request (T728)

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)	7 days	
	Refrigerated	7 days	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
PINW	Q00M0048	Specimen Source	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
PINW	Pinworm Exam, Perianal	Alphanumeric		675-9

 $\ensuremath{\mathsf{LOINC}}\xspace^{\ensuremath{\mathsf{B}}}$ and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87172

Reference Values:

Negative (reported as positive or negative)