

**Overview**
**Method Name**

Varies

**NY State Available**

No

**Specimen**
**Specimen Type**

Varies

**Specimen Required**

Varies

This is a Miscellaneous Referral Test. For specific requirements, please refer to the referred tests list under the catalog tab on the MCL Website (external) or the Referral Catalog (CRM) (internal). If unable to find the correct test, request forms, or specimen requirements, contact Customer Service at 1-800-533-1710 or 6-5700.

NOTE: Provide when ordering

1. Test name
2. Performing lab code
3. Specimen Type
4. Required consent and/or requisitions form including patient specific information, i.e. consent forms, clinical information, family history. Please contact Customer Service at 1-800-533-1710 or 6-5700 for required forms.

**Specimen Minimum Volume**

Varies

**Reject Due To**

Specimens other than	Varies
Anticoagulants other than	NA
Hemolysis	NA
Lipemia	NA
Icteric	NA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

## Clinical and Interpretive

### Reference Values

Test Performed by: Joli Diagnostic

2451 Wehrle Drive

Williamsville, NY 14221

## Performance

### PDF Report

Referral

### Day(s) and Time(s) Test Performed

Varies

### Analytic Time

Varies

### Maximum Laboratory Time

14 days

### Performing Laboratory Location

Joli Diagnostics, Inc.

## Fees and Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### CPT Code Information

Varies

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
ZW173	Misc Joli Diag, Inc.	51991-8

Result ID	Test Result Name	Result LOINC Value
ZT173	Test Name	19145-2
ZR173	Result	19146-0
ZF173	Flag	No LOINC Needed



---

Result ID	Test Result Name	Result LOINC Value
ZV173	Reference Value	19147-8
ZU173	Unit of Measure	No LOINC Needed