

Reporting Title: Heavy Metals, Nails **Performing Location:** Rochester

Necessary Information:

Indicate source of nails (fingernails or toenails), if known

Specimen Requirements:

Supplies: Hair and Nails Collection Kit (T565) Specimen Volume: 0.2 g Collection Instructions:

- 1. Prepare and transport specimen per the instructions in the kit or see Collecting Hair and Nails for Metals Testing.
- 2. Clippings should be taken from all 10 fingernails or toenails.

Specimen Minimum Volume:

0.05 g

Specimen Type	Temperature	Time	Special Container
Nail	Ambient (preferred)		
	Frozen		
	Refrigerated		

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
ASNA	ASNSC	Specimen Source: • Finger • Toe • Finger and Toe	Answer List	Yes
PBNA	PBNSC	Specimen Source: • Finger • Toe • Finger and Toe	Answer List	Yes
HGNA	HGNSC	Specimen Source: • Finger • Toe • Finger and Toe	Answer List	Yes



Heavy Metals, Nails

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
2535	Arsenic, Nails	Numeric	mcg/g	8157-0
	Also used by tests: ASNA			
ASNSC	Specimen Source	nen Source Alphanumeric		31208-2
	Also used by tests: ASNA			
2506	Lead, Nails	Numeric	mcg/g	8202-4
	Also used by tests: PBNA			
PBNSC	Specimen Source	Alphanumeric		31208-2
	Also used by tests: PBNA			
2509	Mercury, Nails	Numeric	mcg/g	8204-0
	Also used by tests: HGNA			
HGNSC	Specimen Source	Alphanumeric		31208-2
	Also used by tests: HGNA			

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
ASNA	Arsenic, Nails			Yes	Yes
PBNA	Lead, Nails			Yes	Yes
HGNA	Mercury, Nails			Yes	Yes

CPT Code Information:

82175- ASNA 83655- PBNA 83825- HGNA



Heavy Metals, Nails

Reference Values:

ARSENIC

0-15 years: Not established > or =16 years: <1.0 mcg/g of nails

LEAD

<4.0 mcg/g of nails Reference values apply to all ages.

MERCURY 0-15 years: Not established > or =16 years: <1.0 mcg/g of nails