

Prothrombin Time (PT), Plasma

# **Reporting Title:** Prothrombin Time (PT), P **Performing Location:** Rochester

#### **Necessary Information:**

Heparin or Coumadin therapy should be noted.

#### **Specimen Requirements:**

Only orderable as part of a profile or reflex. For more information see: ALUPP / Lupus Anticoagulant Profile, Plasma ALBLD / Bleeding Diathesis Profile, Limited, Plasma AATHR / Thrombophilia Profile, Plasma and Whole Blood APROL / Prolonged Clot Time Profile, Plasma ADIC / Disseminated Intravascular Coagulation/Intravascular Coagulation and Fibrinolysis (DIC/ICF) Profile, Plasma

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
PTSEC	Prothrombin Time (PT), P	Numeric	sec	5902-2
INRSC	INR	Numeric		6301-6

LOINC and CPT codes are provided by the performing laboratory.

# Supplemental Report:

No

# **CPT Code Information:**

85610



### **Reference Values:**

Only orderable as part of a profile or reflex. For more information see: ALUPP / Lupus Anticoagulant Profile, Plasma ALBLD / Bleeding Diathesis Profile, Limited, Plasma AATHR / Thrombophilia Profile, Plasma and Whole Blood APROL / Prolonged Clot Time Profile, Plasma ADIC / Disseminated Intravascular Coagulation/Intravascular Coagulation and Fibrinolysis (DIC/ICF) Profile, Plasma

PROTHROMBIN TIME 9.4-12.5 seconds

INTERNATIONAL NORMALIZED RATIO (INR) 0.9-1.1 Standard intensity warfarin therapeutic range: 2.0-3.0 High intensity warfarin therapeutic range: 2.5-3.5