**Test Definition: ECADB** E-Cadherin Immunostain, Bone Marrow, Technical Component

Only

# **Reporting Title:** Bone Marrow E-Cadherin IHC, T Only **Performing Location:** Rochester

MAYO CLINIC LABORATORIES

# **Ordering Guidance:**

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

## **Shipping Instructions:**

Attach the green pathology address label and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

# **Specimen Requirements:**

Supplies: Immunostain Technical Only Envelope (T693) Specimen Type: Tissue Container/Tube: Immunostain Technical Only Envelope Preferred: 2 Unstained positively charged glass slides (25- x 75- x 1-mm) per test ordered; sections 4-microns thick; B5 fixed, decalcified, paraffin-embedded bone marrow tissue.

Acceptable: B5 fixed, decalcified, paraffin-embedded bone marrow tissue block.

#### Forms:

If not ordering electronically, complete, print, and send an Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains Request (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
603212	Bone Marrow E-Cadherin IHC, T Only	Alphanumeric		Bill only; no re

LOINC and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No



Only

## **CPT Code Information:**

88342-TC, primary 88341-TC, if additional IHC

### **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
IHTOI	IHC Initial, Tech Only			No	No
IHTOA IHC Additional, Tech Only				No	No