**Test Definition: SP263** 

BORATORIES Programmed Death-Ligand 1 (PD-L1) (SP263), Semi-Quantitative

Immunohistochemistry, Manual

# **Reporting Title:** PD-L1 (SP263) SemiQuant IHC, Manual **Performing Location:** Rochester

## Ordering Guidance:

MAYO CLINIC

For information on selection of programmed cell death 1-ligand 1 (PD-L1) testing, see https://news.mayocliniclabs.com/pd-I1-by-immunohistochemistry/

### **Shipping Instructions:**

Attach the green pathology address label included in the kit to the outside of the transport container.

### Necessary Information:

A pathology/diagnostic report and a brief history, including primary site of neoplasm, are required.

## **Specimen Requirements:**

Specimen Type: Tissue Supplies: Pathology Packaging Kit (T554) Collection Instructions: Formalin-fixed, paraffin-embedded tissue block; or 3 unstained glass, "positively charged" slides with 4-microns formalin-fixed, paraffin-embedded tissue Additional Information: One slide will be stained with hematoxylin and eosin and returned.

### Forms:

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen: Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains Request (T763) Oncology Test Request (T729)

Specimen Type	Temperature	Time	Special Container
Special	Ambient (preferred)		
	Refrigerated		

# Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
603756	Interpretation	Alphanumeric		59465-5

 MAYOCLINIC
 Test Definition: SP263

 LABORATORIES
 Programmed Death-Ligand 1 (PD-L1) (SP263), Semi-Quantitative

Immunohistochemistry, Manual

Result ID	Reporting Name	Туре	Unit	LOINC®
603757	Participated in the Interpretation	Alphanumeric		No LOINC Needed
603758	Report electronically signed by	Alphanumeric		19139-5
603759	Material Received	Alphanumeric		81178-6
603760	Disclaimer	Alphanumeric		62364-5
603761	Case Number	Alphanumeric		80398-1

LOINC and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

# **CPT Code Information:**

88360