

# **Test Definition: 2SC** S-(2-Succinyl)-Cysteine (2SC) Immunostain, Technical

**Component Only** 

# **Reporting Title:** 2SC IHC, Tech Only **Performing Location:** Rochester

#### Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order, PATHC / Pathology Consultation.

#### **Shipping Instructions:**

Attach the green pathology address label and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

#### **Specimen Requirements:**

Supplies: Immunostain Technical Only Envelope (T693) Specimen Type: Tissue Container/Tube: Immunostain Technical Only Envelope Preferred: 2 unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick. Acceptable: Formalin-fixed, paraffin-embedded (FFPE) tissue block

# Specimen Minimum Volume:

See Specimen Required

#### Forms:

If not ordering electronically, complete, print, and send an Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains Request (T763) with the specimen.

Specimen Type	Temperature Time		Special Container	
TECHONLY	Ambient (preferred)			
	Refrigerated			

# Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
610032	2SC IHC, Tech Only	Alphanumeric		No LOINC Needed

LOINC and CPT codes are provided by the performing laboratory.



# **Supplemental Report:**

No

# **CPT Code Information:**

88342-TC, primary 88341-TC, if additional IHC

# **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
IHTOI	IHC Initial, Tech Only			No	No
ΙΗΤΟΑ	IHTOA IHC Additional, Tech Only			No	No