

GATA-Binding Protein 2, GATA2, Full Gene Analysis, Next-Generation Sequencing, Varies

Reporting Title: GATA2 Gene, Full Gene Analysis

Performing Location: Rochester

#### **Ordering Guidance:**

For cases where the differential diagnosis remains broad, GATA2 may be evaluated as part of a gene panel. See HLHGP / Primary Hemophagocytic Lymphohistiocytosis (HLH) Gene Panel, Varies; SCCNP / Severe Congenital and Cyclic Neutropenia Gene Panel, Varies; or EBLPD / Epstein Barr Virus (EBV) Susceptibility and Lymphoproliferative Disorders Gene Panel, Varies.

Targeted testing for familial variants (also called site-specific or known variants testing) is available for variants identified in the GATA2 gene. See FMTT / Familial Variant, Targeted Testing, Varies. To obtain more information about testing option, call 800-533-1710.

#### **Shipping Instructions:**

Specimen preferred to arrive within 96 hours of collection.

#### **Specimen Requirements:**

Patient Preparation: A previous bone marrow transplant from an allogenic donor will interfere with testing. Call 800-533-1710 for instructions for testing patients who have received a bone marrow transplant.

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Any anticoagulant Specimen Volume: 3 mL Collection Instructions:

1. Invert several times to mix blood.

2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Stability Information: Ambient (preferred) 4 days/Refrigerated

Specimen Type: Skin biopsy

Supplies: Fibroblast Biopsy Transport Media (T115)

Container/Tube: Sterile container with any standard cell culture media (eg, minimal essential media, RPMI 1640). The

solution should be supplemented with 1% penicillin and streptomycin.

Specimen Volume: 4-mm punch

Specimen Stability Information: Refrigerated (preferred)/Ambient

Additional Information: A separate culture charge will be assessed under CULFB / Fibroblast Culture for Biochemical or

Molecular Testing. An additional 3 to 4 weeks is required to culture fibroblasts before genetic testing can occur.

Specimen Type: Cultured fibroblasts

Container/Tube: T-25 flask Specimen Volume: 2 Flasks

Collection Instructions: Submit confluent cultured fibroblast cells from a skin biopsy from another laboratory. Cultured



GATA-Binding Protein 2, GATA2, Full Gene Analysis, Next-Generation Sequencing, Varies

cells from a prenatal specimen will not be accepted.

Specimen Stability Information: Ambient (preferred)/Refrigerated (<24 hours)

Additional Information: A separate culture charge will be assessed under CULFB / Fibroblast Culture for Biochemical or Molecular Testing. An additional 3 to 4 weeks is required to culture fibroblasts before genetic testing can occur.

#### **Specimen Minimum Volume:**

Blood: 1 mL; Skin biopsy or cultured fibroblasts: See Specimen Required

#### Forms:

- 1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available:
- -Informed Consent for Genetic Testing (T576)
- -Informed Consent for Genetic Testing (Spanish) (T826)
- 2. GATA2 Gene Sequencing Patient Information (T811) is recommended.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

## Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
CULFB	CG770	Reason for Referral	Plain Text	No
CULFB	CG899	Specimen	Plain Text	No

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
619803	Test Description	Alphanumeric		62364-5
619804	Specimen	Alphanumeric		31208-2
619805	Source	Alphanumeric		31208-2
619806	Result Summary	Alphanumeric		50397-9
619807	Result	Alphanumeric		82939-0
619808	Interpretation	Alphanumeric		69047-9
619809	Additional Results	Alphanumeric		82939-0



GATA-Binding Protein 2, GATA2, Full Gene Analysis, Next-Generation Sequencing, Varies

Result ID	Reporting Name	Туре	Unit	LOINC®
619810	Resources	Alphanumeric		99622-3
619811	Additional Information	Alphanumeric		48767-8
619812	Method	Alphanumeric		85069-3
619813	Genes Analyzed	Alphanumeric		82939-0
619814	Disclaimer	Alphanumeric		62364-5
619815	Released By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

Supplemental

#### **CPT Code Information:**

81479

88233- Tissue culture, skin, solid tissue biopsy (if appropriate)

88240- Cryopreservation (if appropriate)

#### **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
CULFB	Fibroblast Culture for Genetic Test			No	Yes

## **Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
CULFB	52327	Result Summary	Alphanumeric		50397-9
CULFB	52329	Interpretation	Alphanumeric		69965-2
CULFB	52328	Result	Alphanumeric		82939-0
CULFB	CG770	Reason for Referral	Alphanumeric		42349-1
CULFB	CG899	Specimen	Alphanumeric		31208-2
CULFB	52331	Source	Alphanumeric		31208-2



GATA-Binding Protein 2, GATA2, Full Gene Analysis, Next-Generation Sequencing, Varies

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
CULFB	52332	Method	Alphanumeric		85069-3
CULFB	54625	Additional Information	Alphanumeric		48767-8
CULFB	52333	Released By	Alphanumeric		18771-6

### **Reference Values:**

An interpretive report will be provided