

Hemosiderin, Random, Urine

# **Reporting Title:** Hemosiderin, Random, U **Performing Location:** Rochester

## **Specimen Requirements:**

Container/Tube: Plastic urine container Specimen Volume: 13 mL Collection Instructions: 1. Collect a random urine specimen. 2. No preservative.

# Specimen Minimum Volume:

12 mL

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	7 days	
	Frozen	7 days	
	Ambient	2 hours	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
HSDU2	Hemosiderin, Random, U	Alphanumeric		4644-1

LOINC and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

## **CPT Code Information:**

83070



Hemosiderin, Random, Urine

## **Reference Values:**

Negative