

## **Test Definition: OXYMU**

Oxymorphone Confirmation, Random, Urine

Reporting Title: Oxymorphone Confirmation, U

Performing Location: Rochester

### **Ordering Guidance:**

For situations where chain of custody is required, a Chain-of-Custody Kit (T282) is available. For chain-of-custody testing, order OXYCX / Oxycodone with Metabolite Confirmation, Chain of Custody, Random, Urine.

Additional drug panels and specific requests are available; call 800-533-1710 or 507-266-5700.

### **Additional Testing Requirements:**

If urine creatinine is required or adulteration of the sample is suspected, the following test should also be ordered ADULT / Adulterants Survey, Random, Urine.

### **Specimen Requirements:**

Supplies: Sarstedt 5 mL Aliquot Tube (T914) Collection Container/Tube: Plastic urine container Submission Container/Tube: Plastic, 5-mL tube

Specimen Volume: 3 mL Collection Instructions:

- 1. Collect a random urine specimen.
- 2. No preservative.

Additional Information:

- 1. No specimen substitutions.
- 2. STAT requests are not accepted for this test.
- Submitting <1 mL will compromise our ability to perform all necessary testing.</li>

### Specimen Minimum Volume:

1 mL

### Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	28 days	
	Frozen	28 days	
	Ambient	14 days	



# **Test Definition: OXYMU**

Oxymorphone Confirmation, Random, Urine

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
62622	Oxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	17395-5
42014	Noroxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	90894-7
36024	Oxymorphone Interpretation	Alphanumeric		69050-3

LOINC and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

No

### **CPT Code Information:**

80365 G0480 (if appropriate)

### **Reference Values:**

Negative

Cutoff concentrations by liquid chromatography-tandem mass spectroscopy:

Oxymorphone: 25 ng/mL Noroxymorphone: 25 ng/mL