

## **Test Definition: MSH6**

MSH6 Immunostain, Technical Component Only

Reporting Title: MSH6 IHC, Tech Only

Performing Location: Rochester

### Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

### **Shipping Instructions:**

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

### **Necessary Information:**

If sending normal and tumor blocks; indicate the block number to be stained in performing lab notes (electronic orders) or on the enclosed paperwork (manual orders).

### **Specimen Requirements:**

Supplies: Immunostain Technical Only Envelope (T693)

Specimen Type: Tissue

Supplies: Immunostain Technical Only Envelope (T693) Container/Tube: Immunostains Technical Only Envelope

Preferred: 2 Unstained positively charged glass slides (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.

Acceptable: Formalin-fixed, paraffin-embedded tissue block

#### Forms:

If not ordering electronically, complete, print, and send a Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains Request (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		



# **Test Definition: MSH6**

MSH6 Immunostain, Technical Component Only

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
70825	MSH6 IHC, Tech Only	Alphanumeric		Bill only; no re

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

### **CPT Code Information:**

88342-TC, Primary 88341-TC, If additional IHC

### **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
IHTOI	IHC Initial, Tech Only			No	No
IHTOA	IHC Additional, Tech Only			No	No