

Test Definition: BAP1

BAP1 Immunostain, Technical Component Only

Reporting Title: BAP1 IHC, Tech Only

Performing Location: Rochester

Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

Shipping Instructions:

Attach the green pathology address label and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

Specimen Requirements:

Supplies: Immunostain Technical Only Envelope (T693)

Specimen Type: Tissue

Container/Tube: Immunostain Technical Only Envelope

Preferred: 2 unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.

Acceptable: Formalin-fixed, paraffin-embedded (FFPE) tissue block

Specimen Minimum Volume:

See Specimen Required

Forms:

If not ordering electronically, complete, print, and send a Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains Request (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container		
TECHONLY	Ambient (preferred)				
	Refrigerated				

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
71488	BAP1 IHC, Tech Only	Alphanumeric		Bill only; no re

LOINC and CPT codes are provided by the performing laboratory.



Test Definition: BAP1

BAP1 Immunostain, Technical Component Only

Supp	lemental	Report:
------	----------	---------

No

CPT Code Information:

88342-TC, primary 88341-TC, if additional IHC

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
IHTOI	IHC Initial, Tech Only			No	No
IHTOA	IHC Additional, Tech Only			No	No