

# **Test Definition: FADCF**

Adenosine Deaminase in CSF

Reporting Title: Adenosine Deaminase in CSF

Performing Location: ARUP Laboratories

## **Specimen Requirements:**

Collect CSF in a leak-proof container. Centrifuge specimen at room temperature and send 0.5 mL frozen.

# **Specimen Minimum Volume:**

0.2 mL

Specimen Type	Temperature	Time	Special Container
CSF	Frozen (preferred)	30 days	
	Refrigerated	7 days	
	Ambient	24 hours	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z5892	Adenosine Deaminase, CSF	Alphanumeric		35703-8

LOINC and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

#### **CPT Code Information:**

84311

### **Reference Values:**

0-9 U/L