

**Reporting Title:** Amphotericin B Susceptibility**Performing Location:** U.T. Health Science Ctr San Antonio**Specimen Requirements:**

Container/Tube: Agar slant or other appropriate media

Specimen volume: Organism in pure culture

Complete and submit with specimen:

1. Fungus Testing Lab Request Form

Collection Instructions:

1. Organism must be in pure culture, actively growing.

2. Place specimen in a large infectious container (T146) and label as etiologic agent/infectious substance.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient		

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
FAMBS	Amphotericin B Susceptibility	Alphanumeric		Not Provided

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

87999 - Unlisted Microbiology Procedure (Refer to patient report to apply the appropriate CPT code in place of this unlisted CPT code based on specimen organism.)

87188 - mould-MIC microdilution or agar dilution (if appropriate)

87186 - yeast-MIC microdilution or agar dilution (if appropriate)