MAYO CLINIC LABORATORIES

Test Definition: FBARG

Barley IgG

Overview

Method Name
Enzyme Immunoassay (FEIA)

NY State Available
No

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube:
Preferred: Red top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 0.5 mL Serum
Collection Instructions:
1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Send refrigerate.

Specimen Minimum Volume
See Specimen Required

Reject Due To

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	7 days	
	Frozen	365 days	

Clinical & Interpretive

Clinical Information
Refer to www.eurofins-viracor.com/test-menu

Reference Values
Reference ranges have not been established for food-specific IgG tests.

Interpretation
The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for evaluation by diet eliminations and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of this test.

Performance

PDF Report
No

Day(s) Performed
Monday through Friday

Report Available
5 to 7 days

Performing Laboratory Location
Eurofins Viracor

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information
86001

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FBARG	Barley IgG	51917-3

MAYO CLINIC LABORATORIES

Test Definition: FBARG

Barley IgG

Result ID	Test Result Name	Result LOINC® Value
FBARG	Barley IgG	51917-3