MAYO CLINIC LABORATORIES

Test Definition: FALG

Alpha Lactalbumin IgG

Overview

Method Name
Enzyme Immunoassay (FEIA)

NY State Available
No

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube:
Preferred: Red top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 0.5 mL Serum
Collection Instructions:
1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Send refrigerated.

Specimen Minimum Volume
See Specimen Required

Reject Due To
All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	7 days	
	Frozen	365 days	

Clinical & Interpretive

Clinical Information

Refer to www.eurofins-viracor.com/test-menu/

Reference Values

Reference ranges have not been established for food-specific IgG tests.

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

5 to 7 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86001

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FALG	Alpha Lactalbumin IgG	63418-8

Result ID	Test Result Name	Result LOINC® Value
FALG	Alpha Lactalbumin IgG	63418-8