MAYO CLINIC LABORATORIES

Test Definition: FFTPO

Thrombopoietin (TPO)

Overview

Method Name
Immunoassay

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Specimen Type: Serum
Container/Tube: Red or SST
Specimen Volume: 1 mL
Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum frozen in a plastic vial.

Specimen Minimum Volume
0.5 mL

Reject Due To

Gross Lipemia	Reject
Gross Icterus	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	30 days	

Clinical & Interpretive

Reference Values
7 – 99 pg/mL

Performance

PDF Report
No

Day(s) Performed
Wednesday-first week of every month

Report Available
9 to 37 days

Performing Laboratory Location
Quest Diagnostics Nichols Institute

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was performed using a kit that has not been approved or cleared by the FDA. The analytical performance characteristics of this test have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. This test should not be used for diagnosis without confirmation by other medically established means.

CPT Code Information
83520

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FFTPO	Thrombopoietin (TPO)	43237-7

Result ID	Test Result Name	Result LOINC® Value
FFTPO	Thrombopoietin (TPO)	43237-7