MAYO CLINIC LABORATORIES

Test Definition: FANTI

Anti-IgE

Overview

Method Name
ELISA

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube:
Preferred: Red top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL Serum
Collection Instructions:
1. Centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send refrigerate

Specimen Minimum Volume
Serum: 0.5 mL

Reject Due To

All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Ambient (preferred)	7 days	
	Refrigerated	7 days	
	Frozen	28 days	

Clinical & Interpretive

Clinical Information

Refer to www.eurofins-viracor.com/test-menu

Reference Values
Normal

Performance

PDF Report
No

Day(s) Performed
Wednesday

Report Available
7 to 11 days

Performing Laboratory Location
Eurofins Viracor

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information
83516

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FANTI	Anti-IgE	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FANTI	Anti-IgE	Not Provided