MAYO CLINIC LABORATORIES

Test Definition: FBFPI

Bird Fancier's Precipitin Panel I

Overview

Method Name
Gel Diffusion

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube:
Preferred: Red top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL Serum
Collection Instructions:
1. Centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send refrigerate

Specimen Minimum Volume
Serum: 0.5 mL

Reject Due To

Gross hemolysis	Reject
Gross lipemia	Reject
Gross icterus	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Clinical & Interpretive

Clinical Information
Refer to www.eurofins-viracor.com/test-menu

Reference Values
Negative

Interpretation
The gel diffusion method was used to test this patient's serum for the presence of precipitating antibodies (IgG) to the antigens indicated. These antibodies are serological markers for exposure and immunological sensitization. The clinical significance varies, depending on the history and symptoms.

Performance

PDF Report
No

Day(s) Performed
Tuesday, Friday

Report Available
7 to 11 days

Performing Laboratory Location
Eurofins Viracor

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information
86331 x10

Test Definition: FBFPI

Bird Fancier's Precipitin Panel I

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FBFPI	Bird Fancier's Precipitin Panel I	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z3892	Parakeet Droppings Gel Diffusion	24503-5
Z3894	Parrott Droppings Gel Diffusion	24501-9
Z3890	Cockatiel Droppings Gel Diffusion	23843-6
Z3891	Finch Droppings Gel Diffusion	Not Provided
Z3893	Parakeet Serum Gel Diffusion	24502-7
Z3888	Canary Droppings Gel Diffusion	Not Provided
Z3889	Chicken Serum Gel Diffusion	24498-8
Z3895	Parrott Serum Gel Diffusion	24500-1
Z3897	Pigeon/Dove Serum Gel Diffusion	6733-0
Z3896	Pigeon/Dove Droppings Gel Diffusion	24505-0