MAYO CLINIC LABORATORIES

Test Definition: FUABF

Uric Acid, Body Fluid

Overview

Method Name
Quantitative Spectrophotometry

NY State Available
Yes

Specimen

Specimen Type
Body Fluid

Specimen Required
Specimen Type: Body fluid
Sources: Drain, peritoneal/ascites, pleural or synovial fluid
Container/Tube: Plastic vial
Specimen Volume: 1 mL
Collection Instructions:
1. Centrifuge to remove cellular material.
2. Aliquot 1 mL body fluid into a plastic vial.
3. Send frozen.
4. **Specimen source must be provided**.

Specimen Minimum Volume
0.2 mL

Reject Due To

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Body Fluid	Frozen (preferred)	180 days	
	Ambient	24 hours	
	Refrigerated	5 days	

Clinical & Interpretive

Clinical Information

Refer to https://ltd.aruplab.com/

Reference Values
Refer to report

Performance

PDF Report
No

Day(s) Performed
Monday through Sunday

Report Available
3 to 5 days

Performing Laboratory Location
ARUP Laboratories

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by ARUP Laboratories It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purpose.

CPT Code Information
84560

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FUABF	Uric Acid, Body Fluid	53612-8

Result ID	Test Result Name	Result LOINC® Value
Z4608	Uric Acid, Body Fluid	53612-8
Z4609	Uric Acid Fluid Source	31208-2