MAYO CLINIC LABORATORIES

Test Definition: FCLCN

CLCN1 DNA Sequencing Test

Overview

Method Name
Next Generation Sequencing

NY State Available
Yes

Specimen

Specimen Type
Whole Blood EDTA

Specimen Required
Note: Collection date and informed consent are required.

Specimen Type: Whole blood
Container/Tube: Lavender-top (EDTA)
Specimen Volume: 8 mL
Collection Instructions:
1. Invert several times to mix.
2. Send ambient.

Specimen Minimum Volume
6 mL

Reject Due To

Frozen specimens	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Whole Blood EDTA	Ambient (preferred)	10 days	
	Refrigerated	10 days	

Clinical & Interpretive

Clinical Information

Refer to www.athenadiagnostics.com/view-full-catalog

Reference Values

No mutation detected

A final report will be attached in MayoAccess.

Performance

PDF Report

Referral

Day(s) Performed

Varies

Report Available

21 to 35 days

Performing Laboratory Location

Athena Diagnostics

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information

81406

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FCLCN	CLCN1 DNA Sequencing	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FCLCN	CLCN1 DNA Sequencing	Not Provided